



SIRCAL

APPLICATION BULLETIN

11.1

Why should you put up with the impure Argon and risk your ICP test performance?



When Sircal can help you with this...

Having received a number of requests for a high flowrate version of MP-2000 Rare Gas Purifier to purify Argon for Inductively Coupled Plasma (ICP) tests we are pleased to share with you this good news that Sircal now offers a high flowrate configuration that consists of three MP-2000 connected in parallel to handle flowrate up to 30 Litres per minute. The system setup is as simple as the MP-2000 you are familiar with as all is being taken care of at Sircal factory. Contact us now for any further information or quotation.

Ordering information:

Part No.	Description	Unit
P4084	High Flowrate Rare Gas Purifier, 230V	1 set
P4085	High Flowrate Rare Gas Purifier, 110V	1 set

Contact us:

Sircal Instruments (UK) Limited

Tel: +44 (0) 1342 335 309 Fax: +44 (0) 1342 323 608 Email: sales@sircal.co.uk www.sircal.co.uk

